MERCK 默克

MILLIPLEX®

April 2021 • VOLUME 01

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
# 03

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
P09

Respiratory disease in rhesus macaques inoculated with SARS-CoV-2

P11

Genomic monitoring of SARS-CoV-2 uncovers an Nsp1 deletion variant that modulates type I interferon response


P13

The Gastrointestinal Tract Is an Alternative Route for SARS-CoV-2 Infection in a Nonhuman Primate Model

P15

Lactate Elicits ER-Mitochondrial Mg<sup>2+</sup> Dynamics to Integrate Cellular Metabolism

P17

Adult Neural Stem Cells Are Alerted by Systemic Inflammation through TNF-α Receptor Signalling

P19

Effects of 2-Year Walnut-Supplemented Diet on Inflammatory Biomarkers


P21

Atopic dermatitis microbiomes stratify into ecologic dermatotypes enabling microbial virulence and disease severity

P23

Association between circadian rhythm disruption and polycystic ovary syndrome

P29

Recruited CD68+CD206+ macrophages orchestrate graft immune tolerance to prompt xenogeneic-dentin matrix-based tooth root regeneration

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## IgG Antibodies against SARS-CoV-2 Correlate with Days from Symptom Onset, Viral Load and 1 IL-10

Department of Medicine, University of Virginia Health System, Charlottesville, VA, USA  
 Posted December 7, 2020  
<https://doi.org/10.1101/2020.12.05.20244541>

### Abstract

The emergence of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) has resulted in a pandemic of the respiratory disease coronavirus disease 2019 (COVID-19). Antibody testing is essential to identify persons exposed to the virus and potentially in predicting disease immunity. 183 COVID-19 patients (68 of whom required mechanical ventilation) and 41 controls were tested for plasma IgG, IgA and IgM against the SARS-CoV-2 S1, S2, receptor binding domain (RBD) and N proteins using the MILLIPLEX® SARS-CoV-2 Antigen Panel. Plasma cytokines were concurrently measured using the MILLIPLEX® MAP Human Cytokine/Chemokine/Growth Factor Panel A. As expected, the 183 COVID-19 positive patients had high levels of IgG, IgA and IgM anti-SARS-CoV-2 antibodies against each of the viral proteins. Sensitivity of anti-S1 IgG increased from 60% to 93% one week after symptom onset. S1-IgG and S1-IgA had specificities of 98% compared to the 41 COVID-19 negative patients. The 68 ventilated COVID-19 positive patients had higher antibody levels than the 115 COVID-19 positive patients who were not ventilated. IgG antibody levels against S1 protein had the strongest positive correlation to days from symptom onset. There were no statistically significant differences in IgG, IgA and IgM antibodies against S1 based on age. We found that patients with the highest levels of anti-SARS-CoV-2 antibodies had the lowest viral load in the nasopharynx.

medRxiv  
 THE PREPRINT SERVER FOR HEALTH SCIENCES

CSH Cold Spring Harbor Laboratory  
 BMJ Yale

### IgG Antibodies against SARS-CoV-2 Correlate with Days from Symptom Onset, Viral Load and IL-10


Mary K. Young, Christine Kornmeier, Rebecca M. Carpenter, Nick R. Natale, Jennifer M. Sasson, Michael D. Solga, Amy J. Mathers, Melinda D. Poulter, Xiao Qiang, William A. Petri Jr.  
[doi: https://doi.org/10.1101/2020.12.05.20244541](https://doi.org/10.1101/2020.12.05.20244541)

This article is a preprint and has not been peer-reviewed [what does this mean?]. It reports new medical research that has yet to be evaluated and so should not be used to guide clinical practice.

Abstract Full Text Info/History Metrics Preview PDF

Finally, there was a correlation of high plasma IL-10 with low anti-SARS-CoV-2 antibodies. Anti-SARS-CoV-2 antibody levels, as measured by a novel antigen panel, increased within days after symptom onset, achieving > 90% sensitivity and specificity within one week, and were highest in patients who required mechanical ventilation. Antibody levels were inversely associated with viral load but did not differ as a function of age. The correlation of high IL-10 with low antibody response suggests a potentially suppressive role of this cytokine in the humoral immune response in COVID-19.

Figure 1

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Figure 1: IgG, IgA and IgM antibody response to SARS-CoV-2 S1 increased in ventilated patients.

(a-c) IgG, IgA and IgM antibody responses to SARS-CoV-2 S1 in ventilated COVID-19 positive patients (n=68), not ventilated COVID-19 positive patients (n=115), and COVID-19 negative patients (n=41). \*\*\*\*p<0.0001, \*\*\*p<0.001, \*\*p<0.01, \*p<0.05

Figure 2




Figure 2: IgG, IgA and IgM antibody response to SARS-CoV-2 S1 and age.

(a-c) IgG, IgA and IgM antibody responses to SARS-CoV-2 S1 in patients less than 30 years old (n=12), 30-49 years old (n=53), 50-69 years old (n=70) and greater than 70 years old (n=47).

**Figure 3**

**Figure 4**

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
**Figure 3: Correlation of IgG, IgA and IgM antibodies against SARS-CoV-2 S1 and days from symptom onset.**

(a-c) Correlation of IgG, IgA and IgM antibodies against SARS-CoV-2 S1 and days from symptom onset (168 samples from 123 patients).

**Figure 4: Correlation of IgG Antibodies and SARS-CoV-2 Ct Value.**

(a-d) Correlation of IgG antibodies against SARS-CoV-2 S1, S2, RBD and N and SARS-CoV-2 Ct Value (n=50).

**Figure 5**

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**Figure 5: Correlation of IgG Antibodies and IL-10.**

(a-d) Correlation of anti-S1, S2, RBD and N IgG antibodies and IL-10 in ventilated (black, n=51) and not ventilated (grey, n=40) COVID-19 positive patients.

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**Human Cytokine/Chemokine/ Growth Factor Panel A**

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## Respiratory disease in rhesus macaques inoculated with SARS-CoV-2

COVID-19  
研究

Laboratory of Virology, National Institute of Allergy and Infectious Diseases, National Institutes of Health, Hamilton, MT, USA.

Published: 12 May 2020


<https://doi.org/10.1038/s41586-020-2324-7>

### Abstract

An outbreak of coronavirus disease 2019 (COVID-19), which is caused by a novel coronavirus (named SARS-CoV-2) and has a case fatality rate of approximately 2%, started in Wuhan (China) in December 2019. Following an unprecedented global spread<sup>3</sup>, the World Health Organization declared COVID-19 a pandemic on 11 March 2020. Although data on COVID-19 in humans are emerging at a steady pace, some aspects of the pathogenesis of SARS-CoV-2 can be studied in detail only in animal models, in which repeated sampling and tissue collection is possible. Here we show that SARS-CoV-2 causes a respiratory disease in rhesus macaques that lasts between 8 and 16 days. Pulmonary infiltrates, which are a hallmark of COVID-19 in humans, were visible in lung radiographs. We detected high viral loads in swabs from the nose and throat of all of the macaques, as well as in bronchoalveolar lavages; in one macaque, we observed prolonged rectal shedding. Together, the rhesus macaque recapitulates the moderate disease that has been observed in the majority of human cases of COVID-19. The establishment of the rhesus macaque as a model of COVID-19 will increase our understanding of the pathogenesis of this disease, and aid in the development and testing of medical countermeasures.

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Fig. 3 | Cytokine and chemokine levels in the serum of rhesus macaques infected with SARS-CoV-2. The levels of 23 cytokines and chemokines were determined in serum at different time points after inoculation. Levels are displayed only for those cytokines and chemokines for which statistically significant (one-way analysis of variance) differences were observed compared to levels on the day of inoculation. The lower limit of detection is indicated with a dotted line. Serum samples were analysed in duplicate from each macaque for each time point; n = 8 macaques at 0, 1, and 3 dpi and n = 4 macaques thereafter.

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Non-Human Primate Cytokine/Chemokine Panel (Cat. No. PCYTMG-40K-PX23)

#### Non-Human Primate Cytokine/Chemokine Panel 1

- ☉ (Cat. No. PRCYTOMAG-40K)
- 23 (Cat. No. PCYTMG-40K-PX23)
- 23 (Bulk Cat. No. PRCYMAG40PMX23BK)

sCD40L	IL-12/23 (p40)
G-CSF	IL-13
GM-CSF	IL-15
IFN $\gamma$	IL-17A/CTLA8
IL-1 $\beta$	IL-18
IL-1Ra	MCP-1/CCL2
IL-2	MIP-1 $\alpha$ /CCL3
IL-4	MIP-1 $\beta$ /CCL4
IL-5	TGF $\alpha$
IL-6	TNF $\alpha$
IL-8/CXCL8	VEGF-A
IL-10	

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# Genomic monitoring of SARS-CoV-2 uncovers an Nsp1 deletion variant that modulates type I interferon response



Department of Laboratory Medicine, State Key Laboratory of Biotherapy, West China Hospital, Sichuan University and Collaborative Innovation Center for Biotherapy, Chengdu, China  
 Department of Laboratory Medicine and Department of Pediatric Infectious Diseases, Key Laboratory of Birth Defects and Related Diseases of Women and Children of MOE, West China Second University Hospital, Sichuan University, Chengdu, China  
 Department of Respiratory and Critical Care Medicine, Frontier Science Center of Disease Molecular Network, West China Hospital, Sichuan University, Chengdu, China  
 Department of Clinical Laboratory, Zhongnan Hospital of Wuhan University, Wuhan, Hubei, China


## Abstract

The SARS-CoV-2 virus, the causative agent of COVID-19, is undergoing constant mutation. Here, we utilized an integrative approach combining epidemiology, virus genome sequencing, clinical phenotyping, and experimental validation to locate mutations of clinical importance. We identified 35 recurrent variants, some of which are associated with clinical phenotypes related to severity. One variant, containing a deletion in the Nsp1-coding region ( $\Delta$ 500-532), was found in more than 20% of our sequenced samples and associates with higher RT-PCR cycle thresholds and lower serum IFN- $\beta$  levels of infected patients. Deletion variants in this locus were found in 37 countries worldwide, and viruses isolated from clinical samples or engineered by reverse genetics with related deletions in Nsp1 also induce lower IFN- $\beta$  responses in infected Calu-3 cells. Taken together, our virologic surveillance characterizes recurrent genetic diversity and identified mutations in Nsp1 of biological and clinical importance, which collectively may aid molecular diagnostics and drug design.

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
**Fig 6.** Nsp1 mutants downregulate IFN-1 response (B) Concentration of IFN- $\beta$  in the supernatant of Nsp1-expressing A549 cells. GFP and SARS-CoV Nsp1 were used as negative and positive controls. Each dot represents a single biological replicate. Median and range of 75th percentile are also shown in the boxplots. \*indicates statistical significance between mutant and WT Nsp1, and +indicates significance between GFP controls. Mann-Whitney U test was performed.

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**Figure 7.** Growth kinetics and regulation of IFN-I response in wild-type and Nsp1 deletion virion (G) Concentration of IFN- $\beta$  in the supernatant of infected Calu-3 cells. Mean and SEM are shown. Significance between WT and mutant viruses was indicated; \* $p < 0.05$ , \*\* $p < 0.01$ , and \*\*\* $p < 0.001$ . Multiple t tests were performed.

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### Human Cytokine/Chemokine/ Growth Factor Panel A

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- ☑ (Cat. No. HCYTA-60K)
- 38 (Cat. No. HCYTA-60K-PX38) ♦
- 48 (Cat. No. HCYTA-60K-PX48)


sCD40L	IL-13 ♦
EGF ♦	IL-15 ♦
Eotaxin/CCL11 ♦	IL-17A/CTLA8 ♦
FGF-2/FGF-basic	IL-17E/IL-25 ♦
Fit3 Ligand	IL-17F ♦
Fractalkine/CX3CL1	IL-18 ♦
G-CSF ♦	IL-22 ♦
GM-CSF ♦	IL-27
GRO $\alpha$	IP-10/CXCL10 ♦
IFN $\alpha$ 2 ♦	MCP-1/CCL2 ♦
IFN $\gamma$ ♦	MCP-3/CCL7
IL-1 $\alpha$ ♦	M-CSF ♦
IL-1 $\beta$ ♦	MDC/CCL22
IL-1RA ♦	MIG/CXCL9 ♦
IL-2 ♦	MIP-1 $\alpha$ /CCL3 ♦
IL-3 ♦	MIP-1 $\beta$ /CCL4 ♦
IL-4 ♦	PDGF-AA ♦
IL-5 ♦	PDGF-AB/BB ♦
IL-6 ♦	RANTES/CCL5 ♦▲
IL-7 ♦	TGF $\alpha$
IL-8/CXCL8 ♦	TNF $\alpha$ ♦
IL-9	TNF $\beta$ /Lymphotoxin- $\alpha$ (LTA) ♦
IL-10 ♦	VEGF-A ♦
IL-12 (p40) ♦	
IL-12 (p70) ♦	

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## The Gastrointestinal Tract Is an Alternative Route for SARS-CoV-2 Infection in a Nonhuman Primate Model

COVID-19  
研究

National Kunming High-Level Biosafety Primate Research Center, Institute of Medical Biology, Chinese Academy of Medical Sciences and Peking Union Medical College, Yunnan, China

Available online 9 December 2020  
<https://doi.org/10.1053/j.gastro.2020.12.001>

### Abstract

**BACKGROUND & AIMS:** Gastrointestinal (GI) manifestations have been increasingly reported in patients with coronavirus disease 2019 (COVID-19). However, the roles of the GI tract in severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) infection are not fully understood. We investigated how the GI tract is involved in SARS-CoV-2 infection to elucidate the pathogenesis of COVID-19. **METHODS:** Our previously established nonhuman primate (NHP) model of COVID-19 was modified in this study to test our hypothesis. Rhesus monkeys were infected with an intragastric or intranasal challenge with SARS-CoV-2. Clinical signs were recorded after infection. Viral genomic RNA was quantified by quantitative reverse transcription polymerase chain reaction. Host responses to SARS-CoV-2 infection were evaluated by examining inflammatory cytokines, macrophages, histopathology, and mucin barrier integrity. **RESULTS:** Intranasal inoculation with SARS-CoV-2 led to infections and pathologic changes not only in respiratory tissues but also in digestive tissues. Expectedly, intragastric inoculation with SARS-CoV-2 resulted in the productive infection of digestive tissues and inflammation in both the lung and digestive tissues. Inflammatory cytokines were induced by both types of inoculation with SARS-CoV-2, consistent with the increased expression of CD68. Immunohistochemistry and Alcian blue/periodic acid–Schiff staining showed decreased Ki67, increased cleaved caspase 3, and decreased numbers of mucin-containing goblet cells, suggesting that the inflammation induced by these 2 types of inoculation with SARS-CoV-2 impaired the GI barrier and caused severe infections. **CONCLUSIONS:** Both intranasal and intragastric inoculation with SARS-CoV-2 caused

Colonoscopist Quality to Stratify Surveillance for Adenomas <b>1067</b>	Effect of CRC Screening Strategy on Participation: A Randomized Trial <b>1097</b>	Novel CRC Susceptibility Genes From Transcriptome-Wide Study <b>1164</b>	BMP Signaling Mediated by GREM1 & ISLR: Colorectal Carcinogenesis <b>1224</b>
---	---	--	---

**Gastroenterology**  
 Volume 160 / Number 4 March 2021 [www.gastrojournal.org](http://www.gastrojournal.org)

**Streptococcus thermophilus in Preventing CRC**  
 1179

**also in this issue**


**Cancer Screening During the Coronavirus Disease-2019 Pandemic: A Perspective From NCI's PROSPR Consortium** **999**

**Race, Ethnicity, and Socioeconomic Status Are Associated With Prolonged Time to Treatment After a Diagnosis of Colorectal Cancer: A Large Population-Based Study** **1394**

**aga** American Gastroenterological Association

**IF: 14.88**


pneumonia and GI dysfunction in our rhesus monkey model. Inflammatory cytokines are possible connections for the pathogenesis of SARS-CoV-2 between the respiratory and digestive systems.

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
**Figure 3. Inflammatory responses in the digestive system and virus shedding in rhesus monkeys intranasally challenged with SARS-CoV-2.**

(A) Multiplex assay of inflammatory cytokines in tissue lysates. For each panel, MM-O/N-0, MM-N-1, MM-N-4, MM-N-7, and MM-N-14 are shown in 1 column from left to right. The highest level of cytokine, set as 1, is used to normalize other samples within 1 panel, where the color density represents the relative levels of cytokines.

**Figure 5. Inflammatory responses in the digestive system were induced by SARS-CoV-2 infection.**

(A) Multiplex assay of inflammatory cytokines in digestive tissues. For each panel, MM-O/N-0, MM-O-1, MM-O-4, MM-O-7, and MM-O-14 are shown in 1 column from left to right. The highest level of cytokine, set as 1, is used to normalize the other samples within 1 panel, where the color density represents the relative levels of cytokines.

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
## Lactate Elicits ER-Mitochondrial $Mg^{2+}$ Dynamics to Integrate Cellular Metabolism

Department of Medicine, Center for Precision Medicine, University of Texas Health San Antonio, San Antonio, USA  
Received 11 March 2020, Revised 15 June 2020, Accepted 25 August 2020, Available online 8 October 2020.

Published: October 8, 2020  
<https://doi.org/10.1016/j.cell.2020.08.049>

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
### Abstract

$Mg^{2+}$  is the most abundant divalent cation in metazoans and an essential cofactor for ATP, nucleic acids, and countless metabolic enzymes. To understand how the spatio-temporal dynamics of intracellular  $Mg^{2+}$  ( $iMg^{2+}$ ) are integrated into cellular signaling, we implemented a comprehensive screen to discover regulators of  $iMg^{2+}$  dynamics. Lactate emerged as an activator of rapid release of  $Mg^{2+}$  from endoplasmic reticulum (ER) stores, which facilitates mitochondrial  $Mg^{2+}$  ( $mMg^{2+}$ ) uptake in multiple cell types. We demonstrate that this process is remarkably temperature sensitive and mediated through intracellular but not extracellular signals. The ER-mitochondrial  $Mg^{2+}$  dynamics is selectively stimulated by L-lactate. Further, we show that lactate-mediated  $mMg^{2+}$  entry is facilitated by Mrs2, and point mutations in the intermembrane space loop limits  $mMg^{2+}$  uptake. Intriguingly, suppression of  $mMg^{2+}$  surge alleviates inflammation-induced multi-organ failure. Together, these findings reveal that lactate mobilizes

IF: 38.64

$iMg^{2+}$  and links the  $mMg^{2+}$  transport machinery with major metabolic feedback circuits and mitochondrial bioenergetics.

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## Adult Neural Stem Cells Are Alerted by Systemic Inflammation through TNF- $\alpha$ R Receptor Signalling

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


### Figure 5. TNF- $\alpha$ Receptors Differentially Regulate NSCs In Vitro

(I) Quantification of phosphorylated levels of p38 relative to the total amount of the protein assessed after treatment (1 h) with R1ag or R2ag by Luminex-based multiplex assay. Data are normalized to untreated cells ( $n = 2$ ).

(J) Number of wild-type neurospheres formed in the presence (treatment; left) or after pre-treatment (right) with R2ag and the p38 inhibitor SB203580. Data are represented relative to untreated controls.

All graphs show mean  $\pm$  SEM. \* $p < 0.05$ , \*\* $p < 0.01$ , and \*\*\* $p < 0.001$ .

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### Multi-Pathway (Phosphoprotein) – 9 Plex

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## Effects of 2-Year Walnut-Supplemented Diet on Inflammatory Biomarkers

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Center for Nutrition, Healthy Lifestyle, and Disease Prevention School of Public Health  
Loma Linda University

Published: November 10, 2020  
<https://doi.org/10.1016/j.cell.2020.08.049>

### Abstract

Nut consumption may be associated with lower CVD risk because nuts have a consistent cholesterol-lowering effect. A meta-analysis of 24 randomized controlled trials (RCTs) concluded that, compared with control diets, walnut-enriched diets resulted in significant weighted mean differences in low-density lipoprotein cholesterol (5.5 mg/dl), but had no effect on blood pressure or high-sensitivity C-reactive protein (hs-CRP) (2). Nut-enriched diets also affect endothelial function, with weighted mean differences in flow-mediated dilatation of 0.79% in 8 RCTs (3). These modest salutary effects of nut diets, however, cannot fully account for the lower CVD outcomes observed in prospective studies. Given the prevailing theory that inflammation is a major driver of atherosclerosis, 1 potential mechanism linking nut consumption to reduced CVD might be diminished inflammation.

IF: 20.59

**TABLE 1** Baseline Concentrations and Changes in Circulating Inflammatory Biomarkers by Intervention Group

Biomarker	Walnut Diet (n = 324)	Control Diet (n = 310)	p Value*
GM-CSF, pg/ml			
Baseline	18.2 (16.4 to 20.0)	17.9 (16.0 to 19.9)	0.837
Change	-2.10 (-2.89 to -1.31)	-0.002 (-0.81 to 0.82)	<0.001
Percent change	-11.5 (-15.9 to -7.2)	-0.01 (-4.5 to 4.6)	
IFN- $\gamma$ , pg/ml			
Baseline	16.2 (15.4 to 17.0)	15.2 (14.3 to 16.0)	0.083
Change	-1.35 (-1.81 to -0.89)	-0.09 (-0.55 to 0.38)	<0.001
Percent change	-8.3 (-11.2 to -5.5)	-0.6 (-3.6 to 2.5)	
IL-1 $\beta$ , pg/ml			
Baseline	1.19 (1.12 to 1.26)	1.14 (1.08 to 1.21)	0.354
Change	-0.12 (-0.16 to -0.08)	-0.02 (-0.06 to 0.02)	<0.001
Percent change	-10.1 (-13.4 to -6.7)	-1.8 (-5.3 to 1.8)	
IL-6, pg/ml			
Baseline	2.26 (2.09 to 2.44)	2.16 (1.98 to 2.33)	0.395
Change	-0.19 (-0.30 to -0.08)	-0.01 (-0.12 to 0.10)	0.021
Percent change	-8.4 (-13.3 to -3.5)	-0.5 (-5.6 to 4.6)	
TNF- $\alpha$ , pg/ml			
Baseline	6.09 (5.89 to 6.29)	5.86 (5.64 to 6.08)	0.135
Change	-0.40 (-0.57 to -0.24)	-0.09 (-0.26 to 0.08)	0.009
Percent change	-6.6 (-9.4 to -3.9)	-1.5 (-4.4 to 1.4)	
sE-selectin, ng/ml			
Baseline	49.4 (46.4 to 52.4)	47.1 (44.1 to 50.0)	0.277
Change	-1.66 (-2.67 to -0.65)	0.91 (-0.13 to 1.95)	0.001
Percent change	-3.5 (-5.6 to -1.5)	1.8 (-0.4 to 4.0)	
sICAM-1, ng/ml			
Baseline	1.37 (1.27 to 1.46)	1.24 (1.18 to 1.30)	0.037
Change	0.01 (-0.05 to 0.07)	0.03 (-0.03 to 0.09)	0.674
Percent change	1.5 (-2.9 to 5.1)	3.2 (-1.6 to 8.1)	
sVCAM-1, ng/ml			
Baseline	9.38 (9.14 to 9.63)	9.20 (8.96 to 9.44)	0.302
Change	-0.07 (-0.22 to 0.08)	0.04 (-0.11 to 0.20)	0.305
Percent change	-0.6 (-2.2 to 1.0)	0.7 (-1.1 to 2.3)	
SAA, ng/ml			
Baseline	106.9 (94.1 to 119.7)	106.6 (96.0 to 117.2)	0.970
Change	-2.40 (-11.12 to 6.31)	3.23 (-5.68 to 12.14)	0.377
Percent change	-2.2 (-10.4 to 5.9)	3.0 (-5.3 to 11.4)	
hs-CRP, mg/l			
Baseline	2.91 (2.22 to 3.60)	3.25 (2.62 to 3.88)	0.471
Change	-0.01 (-0.06 to 0.04)	-0.003 (-0.05 to 0.05)	0.903
Percent change	-2.8 (-19.6 to 14.1)	1.9 (-13.2 to 17.2)	

Values are mean (95% confidence interval). \*Obtained by analysis of covariance of the change variables, adjusting for center, age, sex, hypertension, diabetes, education years, smoking, baseline statin use, in-trial statin changes, baseline body mass index, in-trial body mass index change, baseline physical activity, in-trial physical activity changes, and baseline cytokine values.

GM-CSF = granulocyte-monocyte colony stimulating factor; hs-CRP = high-sensitivity C-reactive protein; IFN = interferon; IL = interleukin; SAA = serum amyloid A; sE-selectin = soluble E-selectin; sICAM = soluble intercellular adhesion molecule; sVCAM = soluble vascular cell adhesion molecule; TNF = tumor necrosis factor.

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Table 1 Baseline Concentrations and Changes in Circulating Inflammatory Biomarkers by Intervention Group

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Human CVD Panel 2  
(Cat. No. HCVD2MAG-67K)  
Human High Sensitivity T Cell  
(Cat. No. HSTCMAG-28SK)

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## Atopic dermatitis microbiomes stratify into ecologic dermatotypes enabling microbial virulence and disease severity

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Skin Research Institute of Singapore, Agency of Science Technology and Research Institutes, Singapore


Available online 8 October 2020  
<https://doi.org/10.1016/j.jaci.2020.09.031>

### Abstract


**Background:** Atopic dermatitis (AD) is a common skin disease affecting up to 20% of the global population, with significant clinical heterogeneity and limited information about molecular subtypes and actionable biomarkers. Although alterations in the skin microbiome have been described in subjects with AD during progression to flare state, the prognostic value of baseline microbiome configurations has not been explored. Our aim was to identify microbial signatures on AD skin that are predictive of disease fate.

**Methods:** Nonlesional skin of patients with AD and healthy control subjects were sampled at 2 time points separated by at least 4 weeks. Using whole metagenome analysis of skin microbiomes of patients with AD and control subjects (n = 49 and 189 samples), we identified distinct microbiome configurations (dermatotypes A and B). Blood was collected for immunophenotyping, and skin surface samples were analyzed for correlations with natural moisturizing factors and antimicrobial peptides.

**Results:** Dermatotypes were robust and validated across 2 additional cohorts (63 individuals), with strong enrichment of subjects with AD in dermatotype B. Dermatotype B was characterized by reduced microbial richness, depletion of *Cutibacterium acnes*, *Dermacoccus* and *Methylobacterium* species, individual-specific outlier abundance of *Staphylococcus* species (eg, *S. epidermidis*, *S. capitis*, *S. aureus*), and enrichment in metabolic pathways (eg, branched chain amino acids and arginine biosynthesis) and virulence genes (eg,  $\beta$ -toxin,  $\delta$ -toxin) that defined a pathogenic ecology. Skin surface and circulating host biomarkers exhibited a distinct microbial-associated signature that was further reflected in more severe itching, frequent flares, and increased disease severity in patients harboring the dermatotype B microbiome.

**Conclusion:** We report distinct clusters of microbial profiles that delineate the role of microbiome configurations in AD heterogeneity, highlight a mechanism for ongoing inflammation, and provide prognostic utility toward microbiome-based disease stratification.

**Fig 5.** Systemic and phenotypic differences in dermatotype B patients. A, Clustering based on AD-associated serum biomarker revealed a cluster of dermatotype.

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Human Cytokine/Chemokine Panel I  
(Cat. No. HCYTMAG-60K-PX41)

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## Association between circadian rhythm disruption and polycystic ovary syndrome

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Women's Hospital, School of Medicine, Zhejiang University, Hangzhou, China


Published Jan 04, 2021

<https://doi.org/10.1016/j.fertnstert.2020.10.031>

### Abstract

**Objective:** To explore the association of circadian rhythm disruption with polycystic ovary syndrome (PCOS) and the potential underlying mechanism in ovarian granulosa cells (GCs). **Design:** Multicenter questionnaire-based survey, in vivo and ex vivo studies. **Setting:** Twelve hospitals in China, animal research center, and research laboratory of a women's hospital. **Patients/Animals:** A total of 436 PCOS case subjects and 715 control subjects were recruited for the survey. **In vivo and ex vivo studies** were conducted in PCOS-model rats and on ovarian GCs collected from women with PCOS and control subjects. **Intervention(s):** The PCOS rat model was established with the use of testosterone propionate. **Main Outcome Measure(s):** Assay for transposase-accessible chromatin with high-throughput sequencing (ATAC-seq), RNA sequencing, rhythmicity analysis, functional enrichment analysis. **Result(s):** There was a significant correlation between night shift work and PCOS. PCOS-model rats presented distinct differences in the circadian variation of corticotropin-releasing hormone, adrenocorticotropic hormone, prolactin, and a 4-h phase delay in thyrotropic hormone levels. The motif enrichment analysis of ATAC-seq revealed the absence of clock-related transcription factors in specific peaks of PCOS group, and RNA sequencing ex vivo at various time points over 24 hours demonstrated the differential rhythmic expression patterns of women with PCOS. Kyoto Encyclopedia of Genes and Genomes analysis further highlighted metabolic dysfunction, including both carbohydrate and amino acid metabolism and the tricarboxylic acid cycle. **Conclusion(s):** There is a significant association of night shift work with PCOS, and genome-wide chronodisruption exists in ovarian GCs.

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Figure 1 Circadian variation of peripheral hormone levels in PCOS and control subjects. (E) melatonin, (F) corticotropin-releasing hormone (CRH), (G) adrenocorticotrophic hormone (ACTH), (H) thyrotrophic hormone (TSH), (I) prolactin (PRL), (J) LH/FSH, (K) antimüllerian hormone (AMH). Data presented as mean SEM. Unpaired two tailed Student t test, one-way analysis of variance, and the Bonferroni post test were performed for the data.  $P < .05$  was considered to be statistically significant. PCOS = polycystic ovary syndrome.

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## Successful implantation is associated with a transient increase in serum pro-inflammatory cytokine profile followed by a switch to anti-inflammatory cytokine profile prior to confirmation of pregnancy

The Chinese University of Hong Kong

Published: November 30, 2020


<https://doi.org/10.1016/j.fertnstert.2020.10.031>

### Abstract


**Objective:** To compare the changing peripheral levels of inflammation-related cytokine profile during a 9-day period after blastocyst transfer between women who did and did not conceive. **Design:** Prospective, observational, and longitudinal study. **Setting:** University-affiliated hospital. **Patient(s):** Forty-seven women with infertility who were undergoing single day-5 blastocyst transfer were recruited. **Intervention(s):** This prospective observational and longitudinal study on 47 women with infertility was performed in an in vitro fertilization unit from December 2018 to August 2019. The amounts of a range of cytokines was measured on serial blood samples obtained during a 9-day period after blastocyst transfer. **Main Outcome Measure(s):** Serial blood samples were obtained on the day of embryo transfer, and 3, 6, and 9 days afterward for measurement of serum interferon gamma (IFN-g), tumor necrosis factor alpha, interleukin (IL)-2, IL-4, IL-10, IL-12, IL-13, IL-17, IL-18, and IL-22 using cytometric bead arrays; transforming growth factor beta 1 (TGF-b1) was measured using commercial enzyme-linked immunosorbent assay kits. **Result(s):** The cytokine profile was similar between the women who conceived and those who did not on the day of blastocyst transfer. In women who conceived, IFN-g and IL-17 (pro-inflammatory cytokines) exhibited a transient and significant increase on day 3 after blastocyst transfer, which decreased to the baseline levels by day 6. Meanwhile, IL-10 (anti-inflammatory cytokine) was increased significantly on days 6 and 9, and TGF-b1 (anti-inflammatory cytokine) was increased significantly on day 9 after blastocyst transfer. In women who did not conceive, there was a more pronounced increase in IFN-g and IL-17 (pro-inflammatory cytokines) on day 3, which was sustained on days 6 and 9 without a switch to an anti-inflammatory cytokine profile.

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
Figure 1. A comparison of the changing levels of peripheral cytokine parameters on the day of embryo transfer (ET) and on days 3, 6, and 9 after ET between women who did or did not conceive after ET. The results shown represent mean and standard error for peripheral cytokines and selected ratios on the day of ET (ET t 0) and 3, 6, and 9 days after ET (days ET t 3, t6, and t9, respectively). Y-axis is serum cytokine concentration (pg/mL or ng/mL). NP ?nonpregnant group (n=20); P ?pregnant group (n=27). Unpaired t test was used to compare the pregnant and nonpregnant groups; #P<.05; ##P<.01; ###P<.001. Paired t tests were used to compare the baseline level and sequential changes after ET separately for each group of women; \*P<.05; \*\*P<.01; \*\*\*P<.001; \*\*\*\*P<.0001. IFN= interferon; IL =interleukin; TGF=transforming growth factor; TNF= tumor necrosis factor.

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## Tumor mutation burden in blood predicts benefit from neoadjuvant chemo/ radiotherapy in locally advanced rectal cancer

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
Key laboratory of Carcinogenesis and Translational Research (Ministry of Education), Department of Gastrointestinal Surgery III, Peking University Cancer Hospital & Institute, China

Published: January 2021  
<https://doi.org/10.1016/j.ygeno.2020.10.029>

### Abstract

Distant metastasis has been the major concern of prognosis in patients with locally advanced rectal cancer (LARC). The purpose of this study was to investigate the prognostic value of TMB in blood (bTMB) in LARC patients after receiving neoadjuvant chemoradiotherapy (nCRT) and surgery. Using targeted ctDNA sequencing, we revealed that bTMB level at baseline was positively correlated with recurrence-free survival (RFS). Following nCRT, the patients with decreasing TMB tends to have a longer median RFS. bTMB level after surgery was negatively correlated with RFS. The serum cytokines including IFN $\gamma$ , IFN $\alpha$ 2, IL-1 $\beta$ , IL-2 and MIP-1 $\beta$  were significantly higher in pre-nCRT serum with higher bTMB group than that of lower bTMB group. Clonal evolution analysis showed that the pre- and post-nCRT ctDNAs of most cases had shared mutations. In conclusion, we presume that bTMB could potentially improve pre- and post-treatment risk assessment and facilitate individualized therapy for patients with LARC.

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Fig. 5. The serum cytokine/chemokine levels and bTMB. (A-E) Compare the serum cytokine/chemokine levels between baseline bTMB<10 and bTMB>10 groups. (F-H) Compare the serum cytokine/chemokine levels between post-operation bTMB<10 and bTMB>10 groups.

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## Recruited CD68+CD206+ macrophages orchestrate graft immune tolerance to prompt xenogeneic-dentin matrix-based tooth root regeneration

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

State Key Laboratory of Oral Diseases, West China Hospital of Stomatology, Sichuan University, Chengdu, China

Available online 14 October 2020  
<https://doi.org/10.1016/j.bioactmat.2020.09.029>

### Abstract


Successful regenerative medicine strategies of xenogeneic extracellular matrix need a synergistic balance among inflammation, fibrosis, and remodeling process. Adaptive macrophage subsets have been identified to modulate inflammation and orchestrate the repair of neighboring parenchymal tissues. This study fabricated PPAR $\gamma$ - primed CD68+CD206+ M2 phenotype (M2 $\gamma$ ), and firstly verified their anti-inflammatory and tissue-regenerating roles in xenogeneic bioengineered organ regeneration. Our results showed that Th1-type CD3+CD8+ T cell response to xenogeneic-dentin matrix-based bioengineered root complex (xeno-complex) was significantly inhibited by M2 $\gamma$  macrophage in vitro. PPAR $\gamma$  activation also timely recruited CD68+CD206+ tissue macrophage polarization to xeno-complex in vivo. These subsets alleviated proinflammatory cytokines (TNF- $\alpha$ , IFN- $\gamma$ ) at the inflammation site and decreased CD3+CD8+ T lymphocytes in the periphery system. When translated to an orthotopic nonhuman primate model, PPAR $\gamma$ -primed M2 macrophages immunosuppressed IL-1 $\beta$ , IL-6, TNF- $\alpha$ , MMPs to enable xeno-complex to effectively escape immune-mediated rejection and initiate graft-host synergistic integrity. These collective activities promoted the differentiation of odontoblast-like and periodontal-like cells to guide pulp-dentin and cementum-PDLs-bone regeneration and rescued partially injured odontogenesis such as DSPP and periostin expression. Finally, the regenerated root showed structure-biomechanical and functional equivalency to the native tooth. The timely conversion of M1-to-M2 macrophage mainly orchestrated odontogenesis, fibrogenesis, and osteogenesis, which represents a potential modulator for intact parenchymal-stromal tissue regeneration of targeted organs.

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**Fig. 2. M2 $\gamma$  macrophages suppressed the Th1-type CTL response triggered by xeno-complex in vitro.** (C) Identification of M2 $\gamma$  conditional medium (M2 $\gamma$  medium). The anti-inflammatory and pro-inflammatory cytokines were analyzed using Immunology Multiplex MAP.

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**Fig. 5. PPAR $\gamma$ -primed macrophage inhibited acute inflammation after orthotopic implantation of xeno-complex in Non-human primates.** (C) Th1/Th2-associated cytokines (IL-1 $\beta$ , TNF- $\alpha$ , IL-6, AP-1, IL-4, IL-10). All data are normalized to those of GAPDH mRNA and are presented relative to those of Control, set as 1. All data are mean  $\pm$  SEM ( $n \geq 3$ ), \*\* $p < 0.01$ , \*\*\* $p < 0.001$ , \*\*\*\* $p < 0.0001$ . NS, no significance.

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Non-Human Primate Cytokine/Chemokine Panel 1 (Cat. No. PRCYTOMAG-40K )

### Non-Human Primate Cytokine/Chemokine Panel 1

- (Cat. No. PRCYTOMAG-40K)
- (Cat. No. PCYTMG-40K-PX23)
- (Bulk Cat. No. PRCYMAG40PMX23BK)

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sCD40L	IL-12/23 (p40)
G-CSF	IL-13
GM-CSF	IL-15
IFN $\gamma$	IL-17A/CTLA8
IL-1 $\beta$	IL-18
IL-1Ra	MCP-1/CCL2
IL-2	MIP-1a/CCL3
IL-4	MIP-1 $\beta$ /CCL4
IL-5	TGF $\alpha$
IL-6	TNF $\alpha$
IL-8/CXCL8	VEGF-A
IL-10	

# 03

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28 Jan  
2021

IgG Antibodies Against SARS-CoV-2 Correlate to Days from Symptom Onset in COVID-19 Positive Patients

Available On Demand

### Webinar: IgG Antibodies Against SARS-CoV-2 Correlate to Days from Symptom Onset in COVID-19 Positive Patients

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Mary Young  
Research Technician  
University of Virginia,  
School of Medicine

18 Feb  
2021

### MILLIPLEX® Learning Center is now LIVE!

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18 May  
2021

IL-13 is a Pathogenic Driver of COVID-19

### Webinar: IL-13 is a Pathogenic Driver of COVID-19

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Alexandra Donlan  
University of Virginia

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